


Product Datasheet

Sotrovimab (VIR-7831/GSK4182136) - Research Grade Biosimilar Antibody (orb1238020)

Description	Sotrovimab (VIR-7831/GSK4182136) - Research Grade Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / Spike glycoprotein [SARS-CoV-2]
Target	S
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Human IgG
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
Uniprot ID	P0DTC2
NCBI	P0DTC2
Expiration Date	12 months from date of receipt.



SDS-PAGE Image