

Product Datasheet

SPAR_P03420 Biosimilar Antibody (orb1238032)

Description	SPAR_P03420 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / Respiratory syncytial virus
Target	SPAR_P03420
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Human IgG1
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
Uniprot ID	P03420
NCBI	P03420
Expiration Date	12 months from date of receipt.