

## Product Datasheet

# Anti NS1 of dengue virus - Research Grade Biosimilar (orb1238033)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Anti NS1 of dengue virus - Research Grade Biosimilar   |
| <b>Conjugation</b>     | Unconjugated   |
| <b>Immunogen</b>       | Humanized / Dengue virus   |
| <b>Preservatives</b>   | PBS buffer pH 7.5  |
| <b>Concentration</b>   | batch dependent  |
| <b>Storage</b>         | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term. |
| <b>Note</b>            | For research use only  |
| <b>Isotype</b>         | Human IgG1   |
| <b>Clonality</b>       | Recombinant  |
| <b>Purity</b>          | > 95%  |
| <b>Source</b>          | CHO cells  |
| <b>Uniprot ID</b>      | <b>P14340</b>  |
| <b>NCBI</b>            | <b>P14340</b>  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |