

Product Datasheet

LGR4 Biosimilar Antibody (orb1238034)

Description	LGR4 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / LGR4 (G protein-coupled receptor 48, G protein-coupled receptor GPR48, G-protein coupled receptor 48, leucine-rich repeat containing G protein-coupled receptor 4, Leucine-rich repeat-containing G-protein coupled receptor 4) [Homo sapiens]
Target	LGR4
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Human IgG1
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
Uniprot ID	Q9BXB1
NCBI	Q9BXB1
Expiration Date	12 months from date of receipt.