

Product Datasheet

CD28 Biosimilar Antibody (orb1238035)

Description	CD28 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / CD28 (Tp44) [Homo sapiens]
Target	CD28
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Human IgG4
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
Uniprot ID	P10747
NCBI	P10747
Expiration Date	12 months from date of receipt.