


Product Datasheet

TNFRSF8 Biosimilar Antibody (orb1238037)

Description	TNFRSF8 Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / TNFRSF8 (CD30, D1S166E, Ki-1) [Homo sapiens]
Target	TNFRSF8
Preservatives	PBS buffer pH 7.5
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	Human IgG1
Clonality	Recombinant
Purity	> 95%
Source	CHO cells
Uniprot ID	P28908
NCBI	P28908
Expiration Date	12 months from date of receipt.



SDS-PAGE
Image