


Product Datasheet

BSG Biosimilar Antibody (orb1238041)

Description	BSG Biosimilar Antibody
Conjugation	Unconjugated
Immunogen	Humanized / CD147 [Homo sapiens]
Target	BSG
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG1-kappa
Clonality	Recombinant
Source	CHO cells
Uniprot ID	P35613
NCBI	P35613
Expiration Date	12 months from date of receipt.



SDS-PAGE
Image