
Product Datasheet

AIMF2 Antibody (orb1238077)

Description

AIMF2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC, WB

Immunogen

Anti-AIMF2 antibody (orb1238077) was raised against a peptide corresponding to 14 amino acids near the center of human AIMF2. The immunogen is located within 130-180 amino acids of AIMF2.

Target

AIMF2

Preservatives

AIMF2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

AIMF2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW


Predicted: 41kD Observed: 40 kD

Uniprot ID
[Q9BRQ8](#)
NCBI
[Q9BRQ8](#)
Dilution Range


WB: 2 µg/mL; IHC: 1 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples; Immunofluorescence in human samples and Immunohistochemistry in human, mouse and rat samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



WB Validation in Human Mouse and Rat Cel...



WB Validation in Mouse Skeletal Muscle. ...



WB Validation in Rat Tissues. Loading: 1...