

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody (biotin) (orb1238082)

Description

SARS-CoV-2 (COVID-19)

Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Immunogen	Anti-SARS-CoV-2 (COVID-19) Nucleocapsid antibody (biotin) (orb1238082) was raised against a peptide corresponding to 17 amino acids near the carboxy terminus of SARS-CoV-2 (COVID-19) Nucleocapsid protein. The immunogen is located within 350-400 amino acids of SARS-CoV-2 (COVID-19) Nucleocapsid protein.
Target	N
Preservatives	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody (biotin) is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SARS-CoV-2 COVID-19) Nucleocapsid antibody (biotin) can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
NCBI	QHD43423
Expiration Date	12 months from date of receipt.