

---

## Product Datasheet

### **SARS-CoV-2 (COVID-19) Spike RBD Antibody (biotin) (orb1238083)**

**Description**

SARS-CoV-2 (COVID-19) Spike RBD

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Biotin
<b>Immunogen</b>	Anti-SARS-CoV-2 (COVID-19) Spike RBD antibody (biotin) (orb1238083) was raised against a peptide corresponding to 19 amino acids near the carboxy terminus of SARS-CoV-2 (COVID-19) Spike glycoprotein RBD. The immunogen is located within the last 50 amino acids of SARS-CoV-2 (COVID-19) Spike protein RBD.
<b>Target</b>	S
<b>Preservatives</b>	SARS-CoV-2 (COVID-19) Spike RBD antibody (biotin) is supplied in PBS containing 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	SARS-CoV-2 (COVID-19) Spike RBD antibody (biotin) can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>NCBI</b>	<a href="#">QHD43416</a>
<b>Expiration Date</b>	12 months from date of receipt.