
Product Datasheet

Zika Envelope Antibody (orb1238086)

Description

Zika Envelope Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

Zika Envelope antibody was raised against a peptide corresponding to 14 amino acids near the amino terminus of human Zika Envelope.

Target

Zika Envelope

Preservatives

Zika Envelope Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Zika Envelope antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Zika Envelope antibody can be used for Western blot at 1 - 2 µg/mL.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q32ZE1](#)
NCBI
[AMR39836](#)
Dilution Range

Zika Envelope antibody can be used for Western blot at 1 - 2 µg/mL.

Expiration Date

12 months from date of receipt.