

---

## Product Datasheet

### PUUVsSgp1 Antibody (orb1238091)

**Description**

PUUVsSgp1 Antibody

**Species/Host**

Rabbit

**Reactivity**

Virus

**Conjugation**

Unconjugated

**Tested Applications**

ELISA

**Immunogen**

Puumala virus nucleocapsid antibody was raised against a 19 amino acid peptide near the amino terminus of the Puumala virus nucleocapsid. The immunogen is located within amino acids 130 - 180 of Puumala Virus Nucleocapsid.

**Target**

PUUVsSgp1

**Preservatives**

Puumala virus nucleocapsid antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

Puumala virus nucleocapsid antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

**Note**

For research use only

**Application notes**

Puumala virus nucleocapsid antibody can detect 10ng Puumala virus nucleocapsid peptide in ELISA at 1 µg/ml.

**Isotype**

IgG

**Clonality**

Polyclonal

**Uniprot ID**
[P27313](#)
**NCBI**
[ABO93592](#)
**Dilution Range**

Puumala virus nucleocapsid antibody can detect 10ng Puumala virus nucleocapsid peptide in ELISA at 1 µg/ml.

**Expiration Date**

12 months from date of receipt.