

---

## Product Datasheet

### HTNVsSgp1 Antibody (orb1238093)

**Description**

HTNVsSgp1 Antibody

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Hantaan virus nucleocapsid antibody was raised against an 18 amino acid peptide near the amino terminus of the Hantaan virus nucleocapsid. The immunogen is located within the first 50 amino acids of Hantaan Virus Nucleocapsid.
<b>Target</b>	HTNVsSgp1
<b>Preservatives</b>	Hantaan virus nucleocapsid antibody is supplied in PBS containing 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Hantaan virus nucleocapsid antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
<b>Note</b>	For research use only
<b>Application notes</b>	Hantaan virus nucleocapsid antibody can detect 10ng Hantaan virus nucleocapsid peptide in ELISA at 1 µg/ml.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P05133</a>
<b>NCBI</b>	<a href="#">AFA36178</a>
<b>Dilution Range</b>	Hantaan virus nucleocapsid antibody can detect 10ng Hantaan virus nucleocapsid peptide in ELISA at 1 µg/ml.
<b>Expiration Date</b>	12 months from date of receipt.