
Product Datasheet

HTNVsMgp1 Antibody (orb1238094)

Description

HTNVsMgp1 Antibody

Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Hantaan virus glycoprotein antibody was raised against an 18 amino acid peptide near the amino terminus of the Hantaan virus glycoprotein. The immunogen is located within amino acids 70 - 120 of Hantaan Virus Glycoprotein.
Target	HTNVsMgp1
Preservatives	Hantaan virus glycoprotein antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Hantaan virus glycoprotein antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	Hantaan virus glycoprotein antibody can detect 10ng Hantaan virus glycoprotein peptide in ELISA at 1 µg/ml.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08668
NCBI	BAA05012
Dilution Range	Hantaan virus glycoprotein antibody can detect 10ng Hantaan virus glycoprotein peptide in ELISA at 1 µg/ml.
Expiration Date	12 months from date of receipt.