
Product Datasheet

SNVsSgp1 Antibody (orb1238095)

Description

SNVsSgp1 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

Sin Nombre virus nucleocapsid antibody was raised against a 19 amino acid peptide near the center of the Sin Nombre virus nucleocapsid. The immunogen is located within amino acids 170 - 220 of Sin Nombre Virus Nucleocapsid.

Target

SNVsSgp1

Preservatives

Sin Nombre virus nucleocapsid antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Sin Nombre Nucleocapsid antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

Sin Nombre virus nucleocapsid antibody can detect 10ng Sin Nombre virus nucleocapsid peptide in ELISA at 1 µg/ml.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q9E1J8](#)
NCBI
[AAG03030](#)
Dilution Range

Sin Nombre virus nucleocapsid antibody can detect 10ng Sin Nombre virus nucleocapsid peptide in ELISA at 1 µg/ml.

Expiration Date

12 months from date of receipt.