
Product Datasheet

SNVsMgp1 Antibody (orb1238096)

Description

SNVsMgp1 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

Sin Nombre virus glycoprotein 2 antibody was raised against a 16 amino acid peptide near the center of the Sin Nombre virus glycoprotein. The immunogen is located within amino acids 670 - 720 of Sin Nombre Virus Glycoprotein 2.

Target

SNVsMgp1

Preservatives

Sin Nombre virus glycoprotein 2 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Sin Nombre Glycoprotein 2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

Sin Nombre virus glycoprotein 2 antibody can detect 10ng Sin Nombre virus glycoprotein 2 peptide in ELISA at 1 µg/ml.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q89905](#)
NCBI
[NP_941974](#)
Dilution Range

Sin Nombre virus glycoprotein 2 antibody can detect 10ng Sin Nombre virus glycoprotein 2 peptide in ELISA at 1 µg/ml.

Expiration Date

12 months from date of receipt.