
Product Datasheet

ANDVsMgp1 Antibody (orb1238099)

Description

ANDVsMgp1 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

Andes virus glycoprotein antibody was raised against a 19 amino acid peptide near the amino terminus of the Andes virus glycoprotein. The immunogen is located within amino acids 70 - 120 of Andes Virus Glycoprotein.

Target

ANDVsMgp1

Preservatives

Andes virus glycoprotein antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Andes virus glycoprotein antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

Andes virus glycoprotein antibody can detect 10ng Andes virus glycoprotein peptide in ELISA at 1 µg/ml.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q80DP7](#)
NCBI
[AAO86638](#)
Dilution Range

Andes virus glycoprotein antibody can detect 10ng Andes virus glycoprotein peptide in ELISA at 1 µg/ml.

Expiration Date

12 months from date of receipt.