
Product Datasheet

NP Antibody (orb1238104)

Description

NP Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

RVF virus nucleocapsid antibody was raised against a 17 amino acid peptide near the center of the RVF virus nucleocapsid. The immunogen is located within amino acids 70 - 120 of Rift Valley Fever Virus Nucleocapsid.

Target

NP

Preservatives

Rift Valley Fever Virus Nucleocapsid Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Rift Valley Fever Virus Nucleocapsid antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

RVF virus nucleocapsid antibody can detect 10ng RVF virus nucleocapsid peptide in ELISA at 1 µg/mL.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 27 kDa

Uniprot ID
[Q9QKA6](#)
NCBI
[YP_003848707](#)
Dilution Range

RVF virus nucleocapsid antibody can detect 10ng RVF virus nucleocapsid peptide in ELISA at 1 µg/mL.

Expiration Date

12 months from date of receipt.