
Product Datasheet

B303_sMgp1 Antibody (orb1238106)

Description

B303_sMgp1 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

SFTS Virus HB29 Membrane Glycoprotein antibody was raised against an 18 amino acid synthetic peptide near the center of human SFTS Virus HB29 Membrane Glycoprotein. The immunogen is located within amino acids 570 - 620 of SFTS Virus HB29.

Target

B303_sMgp1

Preservatives

SFTS Virus HB29 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SFTS Virus HB29 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

SFTS Virus HB29 Membrane Glycoprotein antibody can be used for detection of SFTS Virus HB29 Membrane Glycoprotein by ELISA. It will detect 10 ng of free peptide at 1 µg/mL.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[F1BA47](#)
NCBI
[ADZ04471](#)
Dilution Range

SFTS Virus HB29 Membrane Glycoprotein antibody can be used for detection of SFTS Virus HB29 Membrane Glycoprotein by ELISA. It will detect 10 ng of free peptide at 1 µg/mL.