
Product Datasheet

RVFV_sM_gp1 Antibody (orb1238111)

Description

RVFV_sM_gp1 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

RVF virus antibody was raised against a 20 amino acid synthetic peptide near the center of the RVF virus. The immunogen is located within amino acids 360 - 410 of Rift Valley Fever Virus.

Target

RVFV_sM_gp1

Preservatives

Rift Valley Fever Virus Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Rift Valley Fever Virus antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

RVF virus antibody can detect 10ng RVF virus peptide in ELISA at 1 µg/mL.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[P03518](#)
NCBI
[P03518](#)
Dilution Range

RVF virus antibody can detect 10ng RVF virus peptide in ELISA at 1 µg/mL.

Expiration Date

12 months from date of receipt.