
Product Datasheet

Avian Influenza Nonstructural Protein 1 Antibody (orb1238114)

Description

Avian Influenza Nonstructural

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested Applications

ELISA

Immunogen

Avian Influenza Nonstructural Protein 1 antibody was raised against a synthetic peptide corresponding to 14 amino acids at the carboxy terminus of the Avian Influenza Nonstructural Protein 1 protein. Efforts were made to use relatively conserved regions of the viral sequence as the antigen. The immunogen is located within the last 50 amino acids of Avian Influenza Nonstructural Protein 1.

Preservatives

Avian Influenza Nonstructural Protein 1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Avian Influenza Nonstructural Protein 1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Avian Influenza Nonstructural Protein 1 antibody can be used for the detection of the avian influenza nonstructural protein 1 protein from influenza A in ELISA.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q2L6Z3](#)