

Product Datasheet

THAP11 Antibody (orb1238176)

Description	THAP11 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Western blot analysis of THAP11 in human...
Immunogen	THAP11 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human THAP11. The immunogen is located within the last 50 amino acids of THAP11.	
Target	THAP11	
Preservatives	THAP11 antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	THAP11 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.	
Note	For research use only	
Application notes	THAP11 antibody can be used for detection of THAP11 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 35 kDa Observed: 38 kDa	
Uniprot ID	Q96EK4	
NCBI	NP_065190	
Dilution Range	THAP11 antibody can be used for detection of THAP11 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	