
Product Datasheet

TAB3 Antibody (orb1238183)

Description

TAB3 Antibody


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, IHC-P

Immunohistochemistry of TAB3 in human ki...

Immunogen

TAB3 antibody was raised against a 13 amino acid synthetic peptide near the amino terminus of human TAB3. The immunogen is located within amino acids 40 - 90 of TAB3.

Target

TAB3

Preservatives

TAB3 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TAB3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

TAB3 antibody can be used for detection of TAB3 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q8N5C8](#)
NCBI
[NP_690000](#)
Dilution Range

TAB3 antibody can be used for detection of TAB3 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.