


Product Datasheet

StreptII-tag Antibody (orb1238193)

Description	StreptII-tag Antibody
Species/Host	Rabbit
Reactivity	Tag
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	StreptII-tag antibody was raised against the amino acid sequence "gggWSHPQFEK".
Preservatives	StreptII-tag Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	StreptII-tag antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	StreptII-tag antibody can be used for the detection of StreptII by Western blot at 0.25 µg/mL.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 100 kDa Observed: 95 kDa
Dilution Range	StreptII-tag antibody can be used for the detection of StreptII by Western blot at 0.25 µg/mL.
Expiration Date	12 months from date of receipt.



Western blot analysis of 50 ng of StreptII...