
Product Datasheet

SLPI Antibody (orb1238204)

Description	SLPI Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	SLPI antibody was raised against a 17 amino acid synthetic peptide from near the center of mouse SLPI. The immunogen is located within amino acids 60 - 110 of SLPI.
Target	SLPI
Preservatives	SLPI Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SLPI antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	SLPI antibody can be used for detection of SLPI by Western blot at 1-2 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03973
NCBI	P97430
Dilution Range	SLPI antibody can be used for detection of SLPI by Western blot at 1-2 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.

