
Product Datasheet

NA Antibody (orb1238227)

Description

NA Antibody

Species/Host

Mouse

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Mouse monoclonal neuraminidase antibody was raised against a synthetic peptide from the seasonal influenza (H1N1) neuraminidase protein.

Target

NA

Preservatives

Seasonal H1N1 Neuraminidase Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Seasonal H1N1 Neuraminidase monoclonal antibody can be stored at -20°C, stable for one year.

Note

For research use only

Application notes

Neuraminidase antibody can be used for the detection of neuraminidase protein from seasonal influenza A (H1N1) in ELISA. It will also detect 50 ng of free peptide at 1 µg/mL.

Isotype

IgG2b

Clonality

Monoclonal


Uniprot ID
[B1AGS8](#)
NCBI
[ACA33620](#)
Dilution Range

Neuraminidase antibody can be used for the detection of neuraminidase protein from seasonal influenza A (H1N1) in ELISA. It will also detect 50 ng of free peptide at 1 µg/mL.

Expiration Date

12 months from date of receipt.

Seasonal Influenza H1 Specific Antibody



Seasonal influenza A N1 antibody (Cat. N...