

Product Datasheet

DLGAP4 Antibody (orb1238245)

Description	DLGAP4 Antibody	 <p>Western blot analysis of SAPAP4 in SK-N-...</p>
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	SAPAP4 antibody was raised against a 13 amino acid synthetic peptide from near the center of human SAPAP4. The immunogen is located within amino acids 610 - 660 of SAPAP4.	
Target	DLGAP4	
Preservatives	SAPAP4 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	SAPAP4 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	SAPAP4 antibody can be used for detection of SAPAP4 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
Uniprot ID	Q9Y2H0	
NCBI	NP_055717	
Dilution Range	SAPAP4 antibody can be used for detection of SAPAP4 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	