
Product Datasheet

RGPD5 Antibody (orb1238263)

Description

RGPLD5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Immunogen

RGPLD5 antibody was raised against a 23 amino acid synthetic peptide from near the center of human RGPLD5. The immunogen is located within amino acids 1130 - 1180 of RGPLD5.

Target

RGPLD5

Preservatives

RGPLD5 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

RGPLD5 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

RGPLD5 antibody can be used for detection of RGPLD5 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
[Q99666](#)
NCBI
[NP_005045](#)
Dilution Range

RGPLD5 antibody can be used for detection of RGPLD5 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.

 A B


116-

95-

Western blot analysis of RGPLD5 in human ...