

---

## Product Datasheet

### CD274 Antibody (Biotin) (orb1238288)

**Description**

CD274 Antibody (Biotin)

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Biotin

**Tested**

ELISA

**Applications**

A sandwich ELISA was performed using the...

**Immunogen**

PD-L1 antibody was raised against the extracellular domain of human PD-L1.

**Target**

CD274

**Preservatives**

PD-L1 Antibody is supplied in PBS containing 1% rAlbumin and 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

PD-L1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

PD-L1 antibody can be used for as detection antibody in Sandwich Elisa at 1 µg/mL.

**Isotype**

IgG1

**Clonality**

Monoclonal

**Uniprot ID**
[Q9NZQ7](#)
**NCBI**
[NP\\_054862](#)
**Dilution Range**

PD-L1 antibody can be used for as detection antibody in Sandwich Elisa at 1 µg/mL.

**Expiration Date**

12 months from date of receipt.

