
Product Datasheet

PDCL Antibody (orb1238289)

Description	PDCL Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	PDCL antibody was raised against a 16 amino acid peptide near the amino terminus of human PDCL. The immunogen is located within the first 50 amino acids of PDCL.
Target	PDCL
Preservatives	PDCL antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	PDCL antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	PDCL antibody can be used for detection of PDCL by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 33 kDa Observed: 35 kDa
Uniprot ID	Q13371
NCBI	NP_005379
Dilution Range	PDCL antibody can be used for detection of PDCL by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.

