

Product Datasheet

PDCL2 Antibody (orb1238290)

Description

PDCL2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

PDCL2 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human PDCL2. The immunogen is located within the last 50 amino acids of PDCL2.

Target

PDCL2

Preservatives

PDCL2 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PDCL2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

PDCL2 antibody can be used for detection of PDCL2 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 27 kDa Observed: 25 kDa

Uniprot ID

[Q8N4E4](#)

NCBI

[NP_689614](#)

Dilution Range

PDCL2 antibody can be used for detection of PDCL2 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.

