


---

## Product Datasheet

### NLRP14 Antibody (orb1238306)

<b>Description</b>	NLRP14 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	NALP14 antibody was raised against a 16 amino acid synthetic peptide near the amino terminus of human NALP14. The immunogen is located within amino acids 200 - 250 of NALP14.
<b>Target</b>	NLRP14
<b>Preservatives</b>	NALP14 Antibody is supplied in PBS containing 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	NALP14 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
<b>Note</b>	For research use only
<b>Application notes</b>	NALP14 antibody can be used for detection of NALP14 by Western blot at 1 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q86W24</a>
<b>NCBI</b>	<a href="#">NP_789792</a>
<b>Dilution Range</b>	NALP14 antibody can be used for detection of NALP14 by Western blot at 1 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of NALP14 in rat b...