

Product Datasheet

LGI4 Antibody (orb1238330)

Description	LGI4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	LGI4 antibody was raised against a 17 amino acid synthetic peptide from near the carboxy terminus of human LGI4. The immunogen is located within amino acids 470 - 520 of LGI4.
Target	LGI4
Preservatives	LGI4 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	LGI4 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	LGI4 antibody can be used for the detection of LGI4 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N135
NCBI	AAM49552
Dilution Range	LGI4 antibody can be used for the detection of LGI4 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western blot analysis of LGI4 in rat bra...