

Product Datasheet

LGI3 Antibody (orb1238331)

Description	LG13 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	LG13 antibody was raised against a 12 amino acid synthetic peptide from near the center of human LG13. The immunogen is located within amino acids 390 - 440 of LG13.
Target	LG13
Preservatives	LG13 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	LG13 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	LG13 antibody can be used for the detection of LG13 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N145
NCBI	AAQ88483
Dilution Range	LG13 antibody can be used for the detection of LG13 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western blot analysis of LG13 in human b...