

Product Datasheet

JAKMIP1 Antibody (orb1238346)

Description	JAKMIP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Rabbit polyclonal JAKMIP1 antibody was raised against a 16 amino acid peptide near the center of human JAKMIP1. The immunogen is located within amino acids 260 - 310 of JAKMIP1.
Target	JAKMIP1
Preservatives	JAKMIP1 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	JAKMIP1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	JAKMIP1 antibody can be used for detection of JAKMIP1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 69 kDa
Uniprot ID	Q96N16
NCBI	NP_653321
Dilution Range	JAKMIP1 antibody can be used for detection of JAKMIP1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western blot analysis of JAKMIP1 in rat ...