
Product Datasheet

INF2 Antibody (orb1238353)

Description	INF2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	INF2 antibody was raised against a peptide corresponding to 14 amino acids near the amino terminus of human INF2.
Target	INF2
Preservatives	INF2 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	INF2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	INF2 antibody can be used for Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human and rat samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 137 kDa Observed: 160 kDa
Uniprot ID	Q27J81
NCBI	NP_071934
Dilution Range	INF2 antibody can be used for Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human and rat samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.

