


Product Datasheet

GAG Antibody (orb1238365)

Description	GAG Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Rabbit polyclonal HIV-1 p24 antibody was raised against a recombinant full-length HIV-1 p24 protein .
Target	GAG
Preservatives	HIV-1 p24 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	HIV-1 p24 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	P24 antibody can be used for detection of p24 by Western blot at 0.5 - 1 µg/mL.
Isotype	IgG
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	P04591
NCBI	AAB50258
Dilution Range	P24 antibody can be used for detection of p24 by Western blot at 0.5 - 1 µg/mL.
Expiration Date	12 months from date of receipt.



Western blot analysis of p24 using p24 a...