

---

## Product Datasheet

### Hemagglutinin Antibody [4F1F2] (orb1238370)

**Description**

Hemagglutinin Antibody [4F1F2]

**Species/Host**

Mouse

**Reactivity**

Virus

**Conjugation**

Unconjugated

**Tested**

ELISA

**Applications**
**Immunogen**

A peptide corresponding to synthetic peptide corresponding to 14 amino acids in the middle of Hemagglutinin.

**Preservatives**

Hemagglutinin Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

Hemagglutinin monoclonal antibody can be stored at -20°C, stable for one year.

**Note**

For research use only

**Application notes**

Hemagglutinin antibody can be used for the detection of the Hemagglutinin protein from the H5N1 strain of avian influenza A in ELISA. It will detect 10 ng of free peptide at 1 µg/mL. Other applications are pending.

**Isotype**

IgG3

**Clonality**


Monoclonal

**Uniprot ID**
[Q692M2](#)
**NCBI**
[AAT76166](#)
**Dilution Range**

Hemagglutinin antibody can be used for the detection of the Hemagglutinin protein from the H5N1 strain of avian influenza A in ELISA. It will detect 10 ng of free peptide at 1 µg/mL. Other applications are pending.

**Expiration Date**

12 months from date of receipt.



Hemagglutinin antibody at 2 µg/mL speci...