


## Product Datasheet

### Hemagglutinin Antibody [4E11E1] (orb1238371)

<b>Description</b>	Hemagglutinin Antibody [4E11E1]
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant protein corresponding to amino acids 17 - 338 of H5 HA1.
<b>Preservatives</b>	Hemagglutinin Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Hemagglutinin monoclonal antibody can be stored at -20°C, stable for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	Hemagglutinin antibody can be used for detection of H5 HA1 by Western blot at 2.5 µg/mL. Hemagglutinin antibody also specifically recognizes Avian H5N1 influenza virus at 2 µg/mL.
<b>Isotype</b>	IgG1, IgG2b
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q692M2</a>
<b>NCBI</b>	<a href="#">AAT76166</a>
<b>Dilution Range</b>	Hemagglutinin antibody can be used for detection of H5 HA1 by Western blot at 2.5 µg/mL. Hemagglutinin antibody also specifically recognizes Avian H5N1 influenza virus at 2 µg/mL.
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of 25 ng of recomb...