


Product Datasheet

Hemagglutinin Antibody [1E7D8] (orb1238373)

Description	Hemagglutinin Antibody [1E7D8]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A peptide corresponding to 13 amino acids in the middle of the Hemagglutinin protein.
Preservatives	Hemagglutinin Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Hemagglutinin monoclonal antibody can be stored at -20°C, stable for one year.
Note	For research use only
Application notes	Hemagglutinin antibody can be used for the detection of the Hemagglutinin protein from the H5N1 strain of avian influenza A in ELISA. It will detect 10 ng of free peptide at 1 µg/mL. Other applications are pending.
Clonality	Monoclonal
Uniprot ID	Q692M2
NCBI	AAT76166
Dilution Range	Hemagglutinin antibody can be used for the detection of the Hemagglutinin protein from the H5N1 strain of avian influenza A in ELISA. It will detect 10 ng of free peptide at 1 µg/mL. Other applications are pending.
Expiration Date	12 months from date of receipt.



Western blot analysis of (A) 5 ng and (B...