

Product Datasheet

GRTP1 Antibody (orb1238383)

Description	GRTP1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	GRTP1 antibody was raised against a 16 amino acid synthetic peptide from near the center of human GRTP1. The immunogen is located within amino acids 80 - 130 of GRTP1.	Western blot analysis of GRTP1 in SK-N-S...
Target	GRTP1	
Preservatives	GRTP1 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	GRTP1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	GRTP1 antibody can be used for detection of GRTP1 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
Uniprot ID	Q5TC63	
NCBI	EAX09206	
Dilution Range	GRTP1 antibody can be used for detection of GRTP1 by Western blot at 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	