

Product Datasheet

gp120 Antibody (orb1238388)

Description	gp120 Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	gp120 antibody was raised against a peptide corresponding to 15 amino acids near the center of HIV-1 gp120.
Target	gp120
Preservatives	gp120 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	gp120 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	gp120 antibody can be used for Western blot at 1 - 2 µg/mL.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 120 kDa Observed: 120 kDa
Uniprot ID	P04578
NCBI	AAB05604
Dilution Range	gp120 antibody can be used for Western blot at 1 - 2 µg/mL.
Expiration Date	12 months from date of receipt.



Western blot analysis 1 ng of recombinan...