
Product Datasheet

TMC8 Antibody (orb1238409)

Description

TMC8 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Immunogen

EVER2 antibody was raised against a 18 amino acid synthetic peptide from near the amino terminus of human EVER2. The immunogen is located within amino acids 50 - 100 of EVER2.

Target

TMC8

Preservatives

EVER2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

EVER2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

EVER2 antibody can be used for detection of EVER2 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

Uniprot ID
[Q8IU68](#)
NCBI
[AAM44454](#)
Dilution Range

EVER2 antibody can be used for detection of EVER2 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of EVER2 in rat th...