
Product Datasheet

CX3CR1 Antibody (Biotin) (orb1238433)

Description

CX3CR1 Antibody (Biotin)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Biotin

Tested Applications

ELISA

Immunogen

Anti-CX3CR1 antibody (2093) was raised against a peptide corresponding to 20 amino acids near the amino terminus of mature human CX3CR1. The immunogen is located within the first 50 amino acids of CX3CR1.

Target

CX3CR1

Preservatives

CX3CR1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

CX3CR1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 40-44 kD Observed: 50 kD

Uniprot ID
[P49238](#)
NCBI
[NP_001328](#)
Expiration Date

12 months from date of receipt.