


Product Datasheet

ACTB Antibody (Biotin) (orb1238445)

Description	ACTB Antibody (Biotin)
Species/Host	Mouse
Reactivity	Tag
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Biotin-cMyc antibody was raised against a synthetic peptide containing EQKLISEEDL.
Target	ACTB
Preservatives	Biotin-cMyc antibody is supplied in PBS containing 1% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Biotin-cMyc antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	Biotin-cMyc antibody can be used for detection of cMyc by Western blot at 0.5 - 1 µ/ml.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P60709
Dilution Range	Biotin-cMyc antibody can be used for detection of cMyc by Western blot at 0.5 - 1 µ/ml.
Expiration Date	12 months from date of receipt.



Western blot analysis of
(1) 200ng,
(2) ...