
Product Datasheet

HA Antibody (orb1238490)

Description

HA Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA

Applications
Immunogen

Avian Influenza A H7N9 Hemagglutinin antibody was raised against a synthetic peptide corresponding to 11 amino acids near the center of the hemagglutinin [Influenza A virus (A/blue-winged teal/Ohio/566/2006(H7N9))] protein. The immunogen is located within amino acids 260 - 310 of Avian Influenza A H7N9 Hemagglutinin.

Target

HA

Preservatives

Avian Influenza A H7N9 Hemagglutinin Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Avian Influenza A H7N9 Hemagglutinin antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

Avian Influenza A H7N9 Hemagglutinin antibody can be used for detection of Avian Influenza A H7N9 Hemagglutinin by ELISA at 1 µg/ml.

Isotype

IgG

Clonality


Polyclonal

Uniprot ID
[R4NN21](#)
NCBI
[ABS89409](#)
Dilution Range

Avian Influenza A H7N9 Hemagglutinin antibody can be used for detection of Avian Influenza A H7N9 Hemagglutinin by ELISA at 1 µg/ml.

Expiration Date

12 months from date of receipt.



H7N9 Hemagglutinin antibody (Cat. No. or...