

Product Datasheet

HA Antibody (orb1238491)

Description	HA Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Avian Influenza A H7N9 Hemagglutinin antibody was raised against a synthetic peptide corresponding to 11 amino acids near the center of the H7N9 [Influenza A virus (A/Shanghai/02/2013(H7N9))] hemagglutinin protein. The immunogen is located within amino acids 260 - 310 of Avian Influenza A H7N9 Hemagglutinin.
Target	HA
Preservatives	Avian Influenza A H7N9 Hemagglutinin Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Avian Influenza A H7N9 Hemagglutinin antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	Avian Influenza A H7N9 Hemagglutinin antibody can be used for detection of Avian Influenza A H7N9 Hemagglutinin by ELISA at 1 µg/ml.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	R4NN21
NCBI	AGL44438
Dilution Range	Avian Influenza A H7N9 Hemagglutinin antibody can be used for detection of Avian Influenza A H7N9 Hemagglutinin by ELISA at 1 µg/ml.
Expiration Date	12 months from date of receipt.



H7N9 Hemagglutinin antibody (Cat. No. or...