


Product Datasheet

APH1A Antibody (orb1238502)

Description	APH1A Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	APH1 antibody was raised against a 13 amino acid synthetic peptide from near the amino terminus of human APH1a. The immunogen is located within amino acids 40 - 90 of APH1.	Western blot analysis of APH1 in RAW264....
Target	APH1A	
Preservatives	APH1 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	APH1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	APH1 antibody can be used for detection of APH1 by Western blot at 0.5 - 2 µg/mL. Despite its predicted molecular weight, APH1 protein often migrates at aberrant locations in SDS-PAGE. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
Uniprot ID	Q96BI3	
NCBI	AAH08732	
Dilution Range	APH1 antibody can be used for detection of APH1 by Western blot at 0.5 - 2 µg/mL. Despite its predicted molecular weight, APH1 protein often migrates at aberrant locations in SDS-PAGE. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	