
Product Datasheet

TYK2 Antibody (orb1238644)

Description

TYK2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IF, WB

Immunogen

Rabbit polyclonal TYK2 antibody was raised against a 17 amino acid peptide near the amino terminus of human TYK2. The immunogen is located within amino acids 90 - 140 of TYK2.

Target

TYK2

Preservatives

TYK2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TYK2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

TYK2 antibody can be used for detection of TYK2 by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW


Predicted: 131 kDa

Uniprot ID
[P29597](#)
NCBI
[NP_003322](#)
Dilution Range


TYK2 antibody can be used for detection of TYK2 by Western blot at 1 - 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of TYK2 in human s...



Immunofluorescence of TYK2 in A549 cells...