


Product Datasheet

TLX1 Antibody (orb1238749)

Description	TLX1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	TLX1 antibody was raised against a 16 amino acid peptide near the amino terminus of human TLX1. The immunogen is located within the first 50 amino acids of TLX1.
Target	TLX1
Preservatives	TLX1 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	TLX1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	TLX1 antibody can be used for detection of TLX1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in mouse samples and Immunohistochemistry in human samples. All other applications and species not yet tested.
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 36 kDa Observed: 35 kDa
Uniprot ID	P31314
NCBI	NP_005512
Dilution Range	TLX1 antibody can be used for detection of TLX1 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in mouse samples and Immunohistochemistry in human samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western blot analysis of TLX1 in mouse s...



Immunohistochemistry of TLX1 in human sp...